



### Safe, solid, secure: the *Invacare Leo*

The *Invacare Leo* is a scooter designed for individuals who value their independence and wish to get out and about unaided. Safety is a key feature of the **Leo**, which does not detract from its stylish and sporty looks. **Leo** offers individuals the freedom and confidence to enjoy their essential daily outings and leisure excursions. With high standards of design and safety engineering, the **Leo** sets a new market standard.

Trusted Solutions, Passionate People





### For a safe, powerful drive

The *Invacare Leo* offers a smooth, safe drive in a scooter that is easy to control and manoeuvre. The **Leo's** intelligent lighting package includes indicators that automatically switch off and a rear brake light for safe driving. In addition to a visual warning for the battery level display, there is an audible warning (beep) when the battery needs recharging.

### Compact

The **Leo** is a compact mini-scooter that can be easily dismantled for transportation during longer trips or weekend breaks. Easy manoeuvrability for getting around corners in the local shopping centre, museum or leisure complex becomes refreshingly easy. Access to the battery and electronics for maintenance or servicing is much easier with the **Leo** as only one part needs to be removed.

### Reliable

Individuals can trust the **Leo** to take them wherever they need to go. From a trip to the shops or a ride to the park, the **Leo** is sturdy and reliable thanks to its exceptionally high quality electronic engineering and power. These elements are protected from corrosion and

water, (splash protection) which safeguard the vital workings of the scooter against dirt and water. A two-step disengaging lever protects the scooter from free-wheeling if inadvertently knocked out of drive mode. The lighting system is cleverly designed to keep functioning even when one light fails.

## Features and Accessories




#### Finite tiller movement

Individuals can easily adjust the tiller to suit their needs.



#### Brake light

Rear brake lighting warns other people that the scooter is slowing down, thus preventing collisions.



#### Comfort seating

Swivel functionality and sliding rails ensure optimal seating position.



#### Two-step disengaging lever

This unique concept prevents the scooter from accidentally freewheeling.

## Features and Accessories



**Built in anti-splash guards**  
Protects electronics and the transaxle from dirt and water.



**Disassembly into five parts**  
No tools are required to dismantle *Invacare Leo*



**Crutch & cane holder**  
Keeps your crutch or cane stable during the drive.



**Lockable front box**  
Secure box for personal valuables.



**Seatbelt**



**Large mirror**  
Can be mounted left or right.




**Heavy-duty storage cover**  
Keeps the scooter clean and dry.



**Charging off-board**  
Simply remove batteries to charge with the off-board charger kit.









**Seat suspension module**  
Adds extra comfort during driving.



**Windscreen**  
Protects the individuals from the wind.

## Technical data

						
4 wheels	470 mm	410 mm	460 - 560 mm	475 mm	590 mm	1220 mm

						
4 wheels	59 kg (excl. batteries) 83 kg (incl. batteries)	136 kg	2 x 36 Ah	2620 mm	60 mm	70 mm

				
4 wheels	10°	240 W	38 km	8 km/h

\* Seat in maximum fifth rearward locking position from the front.

## Colour of shrouds



Moonstone Silver

Onyx Blue

Ruby Red

## Upholstery



Black vinyl

Invacare Limited  
Pencoed Technology Park -  
Pencoed - Bridgend CF35 5AQ -  
United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS -  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

© 2022 Invacare International Sàrl All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Leo - UK - 04/2022 - 1503875