

# ROBOOTER

R E B O O T Y O U R L I F E

## ROBOOTER X40

TECHNOLOGY MEETS  
FUNCTIONAL ELEGANCE  
A TRUE ORIGINAL



DESIGN  
AWARD  
2022



[www.roobooter.com](http://www.roobooter.com)

# ROBOOTER X40

## BUILT FOR YOU

Ergonomic & Aesthetically Pleasing Design

- Independent shock suspension system
- Memory foam seating
- Auto-folding mechanism
- Smart control with app integration
- Monitoring with fault alerts



Speed  
Settings



Cruise  
Control



Remote  
Drive



Voice  
Settings



Battery  
Monitoring



Gyroscope  
for Stability



## DESIGNED TO GO

Everyday Functionality



Tight turning radius



Omni-wheel built for  
rough terrain



Auto-folds to fit in  
small places



Removable lithium  
battery



Built-in lighting for safety

## PRODUCT PARAMETERS

Endurance Mileage	ABOUT 25KM	Body Size	100.5 X 59 X93 CM
Max Speed	6KM/H	Fold Body Size	82 X 59 X 37 CM
Max Obstacle Height	6CM	Body Net Weight	ABOUT 41.8KG
Ditch Width	10CM	Charging Time	4-6H
Climbing Angle	10°	Front Tire	10-INCH OFF-ROAD OMNIDIRECTIONAL WHEEL
Chassis Net Height	7CM	Rear Tire	10-INCH PNEUMATIC GEAR TIRE
Minimum Rotation Radius	76CM	Side Arms	ANGLE ADJUSTABLE
Seat Height	57.5CM	Maximum Load	120KG
Seat Width	42CM	Battery	TERNARY LITHIUM BATTERY (20AH)
Seat Depth	42CM	Recommended Ambient Temperature	-10°C~50°C
Cushion	42 X 42CM	Waterproof Grade	IPX4