

Please contact Invacare EU Export Customer Services on :-
Tel : +33 247 626980 ; Email : serviceclient_export@invacare.com
 Or visit our website : www.invacare.eu.com



Alber® Erivo®

Folding electric wheelchair with rear wheel drive, mainly for indoor use



Descriptive :





- Programmable **LINX** electronics with reversible, **retractable joystick holder**
- Speed 6 km/h
- 1 x 11 Ah battery as standard, range up to 17 km, 30.15 Ah battery optional or standard on the Erivo Premium. 2 year warranty.
- Black pneumatic rear wheels 12.5", black PU front wheels 8"
- Overall width: 620 mm
- Seat width 450 mm.
- Folding armrests.
- Seat depth 470 mm.
- Backrest height: 450 mm.
- Rear wheel suspension.
- Reflective tape on the backrest.
- Weight of the chair: 27.3 kg without batteries
- Warranty : 2 years except batteries and wear parts

Date :

User name :

Customer Code :

Dealer stamp:

CODE	QUANTITY	DESIGNATION	EUROS	
Alber® Erivo®				
1594813		"Alber® Erivo® - 6,4 km/h , Colour Grey Right-hand joystick, reversible right/left, 11 Ah battery".	2 702,81	
1594814		"Alber® Erivo® Premium - 6,4 km/h , Colour Grey Right hand joystick, reversible right/left, 30,15 Ah battery and road equipment".	3 222,75	
Accessories Alber® Erivo®				
1602792 + 1594974		Dual control with mounting	677,02	
1594975		Road kit (2x front lights, 2x rear lights, 2x white reflectors)	142,18	
1594679		Additional 11 Ah battery (not compatible with the 30.15 Ah battery charger)	639,81	
1594680		Additional battery 30.15 Ah	1000,00	
Think of Bodypoint positioning accessories: belts, harnesses, ankle supports.				
RAMPS				
1654128	<input type="checkbox"/>	Basic ramp 2 telescopic parts in aluminum and non-slip coating. Folded dimensions: 1210 x 194 mm. 8.8 kg per pair - maximum capacity: 270 kg	298,58	 1654128
1654129	<input type="checkbox"/>	Quick2go Compact 3-part aluminum telescopic ramp with non-slip coating Folded dimensions: 870 x 237 mm 13.6 kg per pair - maximum capacity: 260 kg	559,24	 1654129
1654130	<input type="checkbox"/>	Aluminium folding mono ramp Folded dimensions: 1824 x 714 mm 13.9 kg - maximum capacity: 300 kg	317,54	 1654130