

# Invacare® | **Birdie® and Birdie® Compact**



## Transportability and patient comfort go hand-in-hand

With maximum space for the user, both the *Invacare Birdie* and the *Birdie Compact* offer a truly comfortable lift and transfer to or from beds, chairs or even the floor. The lifters are designed to ensure that folding and unfolding can be carried out easily and without the need for tools. In addition, the lifter can be dismantled into two parts if needed – also without tools.



**Yes, you can.®**



**Total user comfort**

Maximizes space for the user even when the hoist is in the highest position. This gives a very comfortable lift and transfer since the user can easily be rotated in the hoist. Consequently, it makes the placing of the user into the bed/wheelchair much easier for the carer. The optimization of space creates a high level of comfort since the user's knees will not be close to the actuator of the hoist. You have, then, a very low risk of injuring the knees of the user.

**Ease of handling: simple to fold, unfold and dismantle without tools**

Handling the lifter has been made as smooth and easy as possible as the hoists can easily be folded and unfolded without the need for tools. Hoists take up minimal space when stored and are easy to push and transport. If required, the mast and boom can be easily detached from the base (no tools required).

**TÜV-certification for maximum security**

For total security, the hoists come with a snap-lock system, thus minimizing the risk of accident. The lifters are CE-certified in line with the Medical Device Directive 93/42/ECC. To ensure

the hoists comply with EN ISO 10535 (the European standard for hoists) the hoist has been rigorously tested both by Invacare and the German TÜV institute, an accredited external body.

**Features and Options**



**Long boom – high user comfort**

The user can easily be rotated and positioned. No risk of injuries to the knees.



**Easy to fold and dismantle**

Ease of transport and storage. No tools required when folding / unfolding or when dismantling the hoist.



**Snap lock**

Optimal security to avoid risk of accidents.



**Emergency stop, emergency lowering and integrated charger**

Manual or electric emergency lowering (compact version features only electric lowering). To charge the hoist, simply plug in.

**Features and Options**



**Click on spreader bar**  
Easy to change the spreader bar.



**Manual leg spread**  
Simply step on the pedal to adjust the leg spread.



**Handle for leg spread**  
Handle can easily be attached to operate the manual leg spread.



**Electrical leg spread**  
Available for *Invacare Birdie* Standard version.



**2 point quick release spreader bars**  
Easy to attach and detach. Available widths: 35 cm, 45 cm, 55 cm.



**Push bar and hand control**  
Superior ergonomics for carers.



**Castors**  
100 mm or 75 mm castors available. (braking castors are shown)



**Scale**  
If needed the scale is available approved, according to the medical standard EN 45501, certified and sealed in order to be used to medicate.



**4 point quick release spreader bars**  
Easy to attach and detach. Available widths: 45 cm and 55 cm.

## Technical data

|                       |        |         |       |                 |       |
|-----------------------|--------|---------|-------|-----------------|-------|
|                       |        |         |       |                 |       |
| <b>Birdie</b>         | 115 cm | 65,5 cm | 42 kg | 46,5 - 168,5 cm | 10 cm |
| <b>Birdie Compact</b> | 100 cm | 53 cm   | 36 kg | 51,5 - 158,5 cm | 10 cm |

|                       |        |        |           |        |        |
|-----------------------|--------|--------|-----------|--------|--------|
|                       |        |        |           |        |        |
| <b>Birdie</b>         | 104 cm | 170 kg | 40 cycles | 2,9 Ah | 140 cm |
| <b>Birdie Compact</b> | 89 cm  | 150 kg | 40 cycles | 2,9 Ah | 107 cm |

All measures are taken with 75 mm castors and 2 point spreaderbar.  
For 100 mm castors add 1,5 cm for height measures and 2 cm for width.

## Frame colour



Light grey /  
Dark grey

Invacare International Sarl  
Route de Cité-Ouest 2 - 1196 Gland  
Switzerland  
Tel: +41 (0)22 354 60 10  
Fax: +41 (0)22 354 60 11  
hqeurope@invacare.com  
www.invacare.eu.com

© 2010 Invacare® International Sarl  
All rights reserved.

All information quoted is believed to be correct at time of print.  
Invacare® reserves the right to alter product specifications without prior consultation.

Birdie - EU - 08/2010

