

Softform Premier Visco



Pressure reducing Visco elastic mattress



The *Invacare* **Softform Premier Visco** is an anti-decubitus mattress featuring a deep visco insert pad, which is high-density and temperature-sensitive to the body, allowing it to conform to the user's shape. This affords exceptional comfort and facilitates pressure reduction by minimising peak pressures.

Pressure ulcer prevention

The temperature-sensitive visco elastic insert pad within the *Invacare* **Softform Premier Visco** warms to conform to the user's shape, maximising the body contact area and thereby facilitating good pressure reduction. The mattress cover is two-way stretch and waterproof, and positively helps to reduce shear and friction forces.

Affords exceptional comfort

The high quality visco elastic foam featured in this mattress delivers an exceptionally comfortable support surface.

Features and Options



Cover is resistant to fluid ingress

Features include two-way stretch waterproof polyurethane cover with high frequency welded seams and a full length flap concealing the zip. Also has a toughened PU coated base to prolong the longevity of the mattress.



Modular Construction

Enables easy parts replacement, increasing product lifetime. Insert pad can be interchanged with regular castellated foam, if required.



Does not require turning

Mattress does not require turning, significantly reducing the manual handling issues.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see www.invacare-eu-export.com/softform-premier-visco

	Width	Length	Height	Total product weight	Max. user weight	Max. cleaning temperature	
							Manufactured to European Standards: BS EN 597-1 Cig, BS EN 597-2 Match, BS7177: 1996 Crib Source 5
Softform Premier Visco	Single size: 660 - 900 mm Double size: 910 - 1370 mm King size: 1380 - 1525 mm	1410 - 2080 mm	152 mm	13 kg	247.6 kg	80°	

Colour

Fuchsia

